

RFTN1 antibody

Catalog Number: 114631

Product name

RFTN1 antibody

Specificity

Human; other species not tested.

Antibody description

RFTN1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Positive WB detected in HeLa cells, Raji cells. Observed molecular weight by Western-blot: 63-70 kDa

Preparation

This antibody was obtained by immunization of RFTN1 recombinant protein (Accession Number: XM_006713069). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:114631(RFTN1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:114631(RFTN1 Antibody) at dilution of 1:50



(under 40x lens)

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114631(RFTN1 Antibody) at dilution of 1:2000