

TNFR1 antibody



Catalog Number: 107631

Product name

TNFR1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TNFR1 Mouse Monoclonal antibody. Positive WB detected in Human brain, HL-60 cells. Positive FC detected in Raji cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of TNFR1 recombinant protein (Accession Number: NM_001065). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107631(TNFR1 antibody) at dilution of 1:1000



1X10⁶ Raji cells were stained with 0.2ug TNFR1 antibody (Catalog No:107631, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.