

# ERLIN1 Polyclonal Antibody



Catalog Number: 164331

## Product name

ERLIN1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ERLIN1

## Preparation

Antigen: A synthetic peptide of human ERLIN1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



Western blot - ERLIN1 Polyclonal Antibody

Western blot analysis of extracts of COLO320 cells, using ERLIN1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.