

IL28B/ IFNL3 antibody



Catalog Number: 111731

Product name

IL28B/ IFNL3 antibody

Specificity

Human; other species not tested.

Antibody description

IL28B/ IFNL3 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of IL28B/ IFNL3 recombinant protein (Accession Number: NM_172139). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111731(IL28B Antibody) at dilution of 1:700