

Anti-KDR antibody



Catalog Number: 142731

Product name

Anti-KDR antibody

Immunogen

[Mouse VEGFR2/Flk-1/CD309/KDR \(His Tag\) recombinant protein](#)

Specificity

Mouse VEGFR2 / Flk-1 / CD309 / KDR

Antibody description

Rabbit Monoclonal to Mouse KDR

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse VEGFR2 / Flk-1 / CD309 / KDR and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



VEGFR2 / Flk-1 / CD309 / KDR Antibody (APC), Rabbit MAb, Flow Cytometry