

RNASEH2C antibody



Catalog Number: 114731

Product name

RNASEH2C antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity


Human; other species not tested.

IHC: 1:20-1:200

Antibody description

RNASEH2C Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human lung tissue, human skin tissue, human spleen tissue. Positive WB detected in human testis tissue, HEK-293 cells, HeLa cells, human spleen tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 22 kDa

Validations




Preparation

This antibody was obtained by immunization of RNASEH2C recombinant protein (Accession Number: NM_032193). Purification method: Antigen affinity purified.

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114731(RNASEH2C antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:114731(RNASEH2C antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:



RNASEH2C antibody

Catalog Number: 114731

Immunohistochemical of paraffin-embedded human brain using Catalog No:114731(RNASEH2C antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-RNASEH2C (IP:Catalog No:114731, 3ug; Detection:Catalog No:114731 1:500) with HEK-293 cells lysate 3200ug.