

# Anti-TMEM161B Rabbit antibody



Catalog Number: 166731

## Product name

Anti-TMEM161B Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse

## Applications

WB, IP

## Antibody description

TMEM161B Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

## Preparation


Antigen: Recombinant protein of human TMEM161B

IP: 1/20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of TMEM161B in K562, HeLa cell lysates using TMEM161B Rabbit pAb(1:1000 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.

## Clonality

Polyclonal