

Exatecan Monoclonal Antibody



Catalog Number: 176731

Product name

Exatecan Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

Exatecan conjugated with OVA. Protein G affinity purified

Formulation

PBS(pH7.4), lyophilized. **Reconstitution: Adding sterile water, prepare a stock solution of 2.0 mg/ml.**

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2b, κ

Applications


ELISA

Dilutions

ELISA: 1:5,000-1:20,000

Validations


Mouse anti-Exatecan Antibody ELISA
0.1 μ g/ml



ELISA Binding Assay of Anti-Exatecan Antibody to different proteins.

ELISA analysis of Anti-Exatecan Monoclonal Antibody with Exatecan-OVA, Exatecan-BSA, OVA and BSA. Anti-Exatecan Monoclonal Antibody can bind to Exatecan-BSA or Exatecan-OVA but not carrier proteins.

anti-Exatecan mouse mAb MM23L



ELISA Binding Assay of Anti-Exatecan Antibody

Immobilized Exatecan ADC at 2 μ g/mL (100 μ L/well) can bind anti-Exatecan mouse monoclonal antibody.