

# Human Carboxypeptidase E/CPE (Fc Tag) recombinant protein



Catalog Number: 500041

## General Information

### Gene Name Synonym

Carboxypeptidase H; Enkephalin convertase;  
Prohormone-processing carboxypeptidase

### Protein Construction

A DNA sequence encoding the human carboxypeptidase E (CPE) precursor (NP\_001864.1) (Met 1-Ser 453) was expressed with the C-terminal fused Fc region of human IgG1.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu 26

### Molecular Mass

The recombinant human CPE/Fc is a disulfide-

linked homodimeric protein. The reduced monomer consists of 666 amino acids and predicts a molecular mass of 74.6 kDa. As a result of glycosylation, the apparent molecular mass of rhCPE/Fc monomer is approximately 85-90 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 100mM Glycine, 10mM NaCl, 50mM Tris, pH 7.5

- 5 % trehalose and mannitol are added as protectants before lyophilization.
- Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human Carboxypeptidase E / CPE Protein (Fc Tag) SDS-PAGE