

# Cynomolgus Betacellulin / BTC (Fc Tag) recombinant protein



Catalog Number: 502041

## General Information

### Protein Construction

A DNA sequence encoding the cynomolgus BTC (XP\_001101880.1) (Met1-Tyr111) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Cynomolgus

### Expression Host

Human Cells

## QC Testing

### Activity

Measured in a cell proliferation assay using Balb/3T3 mouse embryonic fibroblast cells. The ED<sub>50</sub> for this effect is typically 5-40 ng/mL.

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Asp 32

### Molecular Mass

The recombinant cynomolgus BTC is a disulfide-linked homodimer. The reduced monomer comprises 321 amino acids and has a calculated molecular mass of 36 KDa. The apparent molecular mass of the protein is approximately 46 KDa respectively in SDS-PAGE.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Cynomolgus BTC / Betacellulin Protein (Fc Tag)  
SDS-PAGE