

Human METAP2/MNPEP/P67EIF2 (His tag) recombinant protein



Catalog Number: 518041

General Information

Gene Name Synonym

Methionine aminopeptidase 2; MAP 2; MetAP 2; p67; p67eIF2; Peptidase M; METAP2; MAP2

Protein Construction

A DNA sequence encoding the Human METAP2/MNPEP/P67EIF2 Ala2-Tyr478 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human METAP2/MNPEP/P67EIF2 (His tag)
recombinant protein