

# Human GALNT7 (His tag) recombinant protein



Catalog Number: 519041

## General Information

### Gene Name Synonym

PolypeptideGalNAc-transferase7,Protein-UDP acetylgalactosaminyltransferase7,UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 7, GALNT7.

### Protein Construction

A DNA sequence encoding the Human GALNT7 Pro30-Val657 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50mM TrisHCl, 10mM reduced Glutathione, PH8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human GALNT7 (His tag) recombinant protein