

Mouse MIF (His Tag) recombinant protein



Catalog Number: 503441

General Information

Gene Name Synonym

Delayed early response protein 6; Glycosylation-inhibiting factor; L-dopachrome isomerase; L-dopachrome tautomerase; Phenylpyruvate tautomerase

Protein Construction

A DNA sequence encoding the mouse MIF (P34884) precursor (Pro 2-Ala 115) was fused with a polyhistidine tag at the C-terminus and a signal peptide at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Pro 2

Molecular Mass

The secreted recombinant mouse MIF consists of 125 amino acids and has a calculated molecular mass of 13.8 kDa. As a result of glycosylation, the apparent molecular mass of the recombinant protein is approximately 18 and 22 kDa in SDS-PAGE under reducing conditions due to different glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse MIF / Migration Inhibitory Factor Protein (His Tag) SDS-PAGE