

Anti-LILRB5 antibody

Catalog Number: 176441

Product name

Anti-LILRB5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to LILRB5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the n-terminal region of human lilrb5. aa range:21-70

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of HepG2 cells using CD85c Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2.; Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:100