

# Anti-TRAF6BP/TAX1BP1 Rabbit antibody



Catalog Number: 167441

## Product name

Anti-TRAF6BP/TAX1BP1 Rabbit antibody

WB: 1/1000

IHC: 1/20

## Specificity

Human, Mouse, Rat

ICC/IF: 1/20

## Antibody description

TRAF6BP/TAX1BP1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human TRAF6BP

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, ICC/IF

## Dilutions

## Validations



Western blot detection of TRAF6BP/TAX1BP1 in K562, HeLa cell lysates using TRAF6BP/TAX1BP1 Rabbit pAb(1:1000 diluted). Predicted band size: 91kDa. Observed band size: 91kDa.



Immunohistochemistry of TRAF6BP in paraffin-embedded Human lung cancer tissue using TRAF6BP Rabbit pAb at dilution 1/20