

# Collagen Type XVI antibody



Catalog Number: 109441

## Product name

Collagen Type XVI antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Collagen Type XVI Rabbit Polyclonal antibody.  
Positive WB detected in human placenta tissue.  
Observed molecular weight by Western-blot:  
158kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001856).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109441(COL16A1-Specific antibody) at dilution of 1:1500