

Anti-NGFR/P75 antibody

Catalog Number: 100841

Product name

Anti-NGFR/P75 antibody

Specificity

Human NGFR / P75 / CD271

Antibody description

Rabbit polyclonal to NGFR/P75

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human NGFR / P75 / CD271, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 10-20 µg/ml


IP: 0.5-2 uL/mg of lysate

Validations



Items	Lanes	A	B
Sample (whole cell lysate)		K562	Jurkat
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	10-20 µg/ml		
Secondary Antibody	Dylight 800 Labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 75 kDa Observed band size : 75 kDa		

NGFR / P75 / CD271 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



Items	Lanes	A
Sample (whole cell lysate)		Jurkat
Sample quantity		0.5 mg
IP antibody quantity		2 µg
Immunomagnetic beads Protein G		60 µg
Gel	13% SDS-PAGE reducing gel	
Primary antibody	K1H-SMCC-sinoA1003 antibody at 10 µg/ml	
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution	

NGFR / P75 / CD271 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation