

# Galectin-7 antibody



Catalog Number: 110841

## Product name

Galectin-7 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Galectin-7 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue, mouse small intestine tissue. Positive IHC detected in human cervix tissue. Observed molecular weight by Western-blot: 14-15kd

## Preparation

This antibody was obtained by immunization of Galectin-7 recombinant protein (Accession Number: NM\_001042507). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110841(Galectin-7 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110841(Galectin-7 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110841(Galectin-7

# Galectin-7 antibody



Catalog Number: 110841

---

antibody) at dilution of 1:100 (under 40x lens)