

SUPT16H antibody



Catalog Number: 115841

Product name

SUPT16H antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

SUPT16H Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007192). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115841(SUPT16H antibody) at dilution of 1:300




Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:115841(SUPT16H antibody) at dilution of 1:100 (under 10x lens)



SUPT16H antibody

Catalog Number: 115841

Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:115841(SUPT16H antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-SUPT16H (IP:Catalog No:115841, 4ug; Detection:Catalog No:115841 1:1000) with mouse brain tissue lysate 3440ug.