

# Anti-phospho-Rb (Thr356) Rabbit antibody



Catalog Number: 167841

## Product name

Anti-phospho-Rb (Thr356) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P, IP

## Antibody description

phospho-Rb (Thr356) Rabbit monoclonal antibody

## Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr356 of human Rb

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Rb (phospho-Thr356) in HL-60,U2OS,C6 using Rb (phospho-Thr356) Rabbit mAb(1:1000 diluted)