

Mouse IFNGR2 (His Tag) recombinant protein



Catalog Number: 501141

General Information

Gene Name Synonym

Interferon gamma receptor 2; Protein Ifngr2

Protein Construction

A DNA sequence encoding the extracellular domain of mouse IFNGR2 (NP_032364.1) (Met 1-Val 243) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 20

Molecular Mass

The secreted recombinant mouse IFNGR2 comprises 234 amino acids and has a calculated molecular mass of 26.7 kDa. As a result of glycosylation, the apparent molecular mass of the recombinant protein is approximately 40-45 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IFNGR2 Protein (His Tag) SDS-PAGE