

RGS5 Polyclonal Antibody



Catalog Number: 165141

Product name

RGS5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RGS5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human RGS5 (NP_003608.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - RGS5 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using RGS5 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - RGS5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using RGS5 antibody at dilution of 1:100 (40x



lens).

Immunohistochemistry - RGS5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using RGS5 antibody at dilution of 1:100 (40x lens).