

CD7 antibody

Catalog Number: 107141

Product name

CD7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD7 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue.

Preparation

This antibody was obtained by immunization of CD7 recombinant protein (Accession Number: NM_006137). Purification method: .

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:100-1:400

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107141 (CD7 Antibody) at dilution of 1:200 (under 10x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107141 (CD7 Antibody) at dilution of 1:200 (under 40x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).