

# Creatine Kinase antibody



Catalog Number: 109541

---

## Product name

Creatine Kinase antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Creatine Kinase Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 43kd

## Preparation

This antibody was obtained by immunization of Creatine Kinase recombinant protein (Accession Number: NM\_001824). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109541 (CKM antibody) at dilution of 1:1000