

plasminogen antibody



Catalog Number: 113941

Product name

plasminogen antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

plasminogen Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 97kd

Preparation

This antibody was obtained by immunization of plasminogen recombinant protein (Accession Number: NM_000301). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113941(PLG antibody) at dilution of 1:1500