

Anti-CYP4X1 antibody



Catalog Number: 175941

Product name

Anti-CYP4X1 antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to CYP4X1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cytochrome p450 4x1. aa range:251-300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using CYP4X1 Polyclonal Antibody