

ZMYND10 antibody

Catalog Number: 116941

Product name

ZMYND10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZMYND10 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IP detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of ZMYND10 recombinant protein (Accession Number: NM_015896). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116941(ZMYND10 antibody) at dilution of



1:400

IP Result of anti-ZMYND10 (IP:Catalog No:116941, 4ug; Detection:Catalog No:116941 1:500) with SH-SY5Y cells lysate 1600ug.