

TOMM5 antibody



Catalog Number: 116241

Product name

TOMM5 antibody

Specificity

Human; other species not tested.

Antibody description

TOMM5 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HepG2 cells. Observed molecular weight by Western-blot: 6-10 kDa

Preparation

This antibody was obtained by immunization of TOMM5 recombinant protein (Accession Number: NM_001001790). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116241(TOMM5 Antibody) at dilution of 1:600