

DGCR8 antibody



Catalog Number: 107241

Product name

DGCR8 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DGCR8 Mouse Monoclonal antibody. Positive WB detected in K562 cells, BxPC-3 cells, HepG2 cells, human brain tissue, K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of DGCR8 recombinant protein (Accession Number: BC009323). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107241(DGCR8 antibody) at dilution of 1:1000



IP Result of anti-DGCR8 (IP:Catalog No:107241, 5ug; Detection:Catalog No:107241 1:500) with K-562 cells lysate 3440ug.