

FGFR1OP2 antibody



Catalog Number: 110641

Product name

FGFR1OP2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FGFR1OP2 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of FGFR1OP2 recombinant protein (Accession Number: NM_001171888). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110641(FGFR1OP2 antibody) at dilution of 1:500