

# RGS11 antibody



Catalog Number: 114641

## Product name

RGS11 antibody

## Specificity

Human, Rat, Mouse; other species not tested.

## Antibody description

RGS11 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, rat brain tissue. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of RGS11 recombinant protein (Accession Number: NM\_183337). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114641(RGS11 Antibody) at dilution of 1:1000