

Anti-LTA antibody



Catalog Number: 176641

Product name

Anti-LTA antibody

Ig Type

Mouse IgG1

Specificity

Human, Rat

Applications

WB

Antibody description

Mouse monoclonal antibody to LTA

Dilutions

ELISA: 1:10,000 WB: 1:500

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 25% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Validations




Fig1: Western blot analysis on recombinant LT-alpha (TNF- β) using anti- LT-alpha monoclonal antibody.