

Anti-DGKA Rabbit antibody



Catalog Number: 167641

Product name

Anti-DGKA Rabbit antibody

Specificity

Human

Antibody description

DGKA Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human DGKA

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of DGKA in HeLa,A549 cell lysates using DGKA Rabbit mAb(1:1000 diluted).Predicted band size:83kDa.Observed band size:83kDa.