

# PUS1 antibody



Catalog Number: 114341

## Product name

PUS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PUS1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HEK-293 cells, HeLa cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of PUS1 recombinant protein (Accession Number: NM\_025215). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114341(PUS1 Antibody) at dilution of 1:600



Immunofluorescent analysis of MCF-7 cells using Catalog No:114341(PUS1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)