

# Human RELT/TNFRSF19L (Fc tag) recombinant protein



Catalog Number: 519341

## General Information

### Gene Name Synonym

Tumor necrosis factor receptor superfamily member 19L, TNFRSF19L, Receptor expressed in lymphoid tissues, RELT

### Protein Construction

A DNA sequence encoding the Human RELT/TNFRSF19L Ser26-Ala160 is expressed with a Fc tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.2.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human RELT/TNFRSF19L (Fc tag) recombinant protein