

# IL36RN antibody

Catalog Number: 111741

## Product name

IL36RN antibody

## Specificity

Human; other species not tested.

## Antibody description

IL36RN Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of IL36RN recombinant protein (Accession Number: NM\_012275). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111741(IL1F5 Antibody) at dilution of 1:1000