

ND2 antibody

Catalog Number: 113051

Product name

ND2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ND2 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in mouse brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 42-44 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: GQ398480). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113051(ND2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:113051(ND2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:113051(ND2 antibody) at dilution of 1:50 (under 40x lens)