

# POLR3F antibody



Catalog Number: 114051

---

## Product name

POLR3F antibody

## Immunogen

[Human POLR3F Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

POLR3F Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, BxPC-3 cells, HeLa cells, mouse skin tissue. Observed molecular weight by Western-blot: 36-39kd

## Preparation

This antibody was obtained by immunization of POLR3F recombinant protein (Accession Number: NM\_006466). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114051(POLR3F antibody) at dilution of 1:400