

# Anti-FABP1 Rabbit antibody



Catalog Number: 167051

## Product name

Anti-FABP1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FABP1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human FABP-1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20-1/100

IP: 1/20

## Validations



Western blot detection of FABP-1 in K562, Rat Brain, 3T3 cell lysates using FABP-1 Rabbit pAb(1:1000 diluted). Predicted band size: 48kDa. Observed band size: 45kDa.