

Anti-p16 ARC Rabbit antibody



Catalog Number: 169051

Product name

Anti-p16 ARC Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

p16 ARC Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human p16 ARC

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

1:2000

1:20

1:20

Validations



Western blot detection of p16 ARC in Rat Brain, C6, 3T3, HeLa cell lysates using p16 ARC Rabbit mAb(1:1000 diluted). Predicted band size: 16kDa. Observed band size: 16kDa.