

UBE2O Polyclonal Antibody



Catalog Number: 163451

Product name

UBE2O Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to UBE2O

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1073-1292 of human UBE2O (NP_071349.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations




Immunohistochemistry - UBE2O Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using UBE2O antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - UBE2O Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using UBE2O antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - UBE2O Polyclonal

UBE2O Polyclonal Antibody



Catalog Number: 163451

Antibody

human stomach using UBE2O antibody at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded