

Anti-MLH1 (2B3) Mouse antibody



Catalog Number: 168451

Product name

Anti-MLH1 (2B3) Mouse antibody

Specificity

Human

Antibody description

MLH1 (2B3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human MLH1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/500

Validations



Western blot detection of MLH1 in HeLa,293T and Jurkat cell lysates using MLH1 mouse mAb (1:500 diluted).Predicted band size:85KDa.Observed band size:85KDa.