

# Anti-MLH1 (2B3) Mouse antibody



Catalog Number: 168451

## Product name

Anti-MLH1 (2B3) Mouse antibody

## Ig Type

Mouse IgG2b

## Specificity

Human

## Applications

WB

## Antibody description

MLH1 (2B3) Mouse monoclonal antibody

## Dilutions

WB: 1/500

## Preparation

Antigen: Purified recombinant human MLH1 protein fragments expressed in E.coli.

## Validations

## Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of MLH1 in HeLa, 293T and Jurkat cell lysates using MLH1 mouse mAb (1:500 diluted). Predicted band size: 85KDa. Observed band size: 85KDa.