

Mouse Frizzled-1/FZD1 (His Tag) recombinant protein



Catalog Number: 500851

General Information

Protein Construction

A DNA sequence encoding the mouse FZD1 (NP_067432.2) extracellular domain (Met 1-His 248) was expressed with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Val 69

Molecular Mass

The secreted recombinant mouse FZD1 pro form consists of 191 amino acids and has a predicted molecular mass of 21 kDa. In SDS-PAGE under

reducing conditions, the apparent molecular mass of rm FZD1 is approximately 35-40 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Frizzled-1 / FZD1 Protein (His Tag) SDS-PAGE