

GALNTL5 antibody



Catalog Number: 110851

Product name

GALNTL5 antibody

Specificity

Human; other species not tested.

Antibody description

GALNTL5 Rabbit Polyclonal antibody. Positive IHC detected in human lung tissue. Positive WB detected in human placenta tissue, human testis tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of GALNTL5 recombinant protein (Accession Number: XM_006715862). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110851(GALNTL5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung using Catalog No:110851(GALNTL5 antibody) at dilution of 1:100 (under 40x lens)