

INSR antibody

Catalog Number: 111851

Product name

INSR antibody

Specificity

Human; other species not tested.

Antibody description

INSR Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HEK-293 cells, HeLa cells, HT-1080 cells, MCF7 cells. Observed molecular weight by Western-blot: 190 kDa, 135 kDa, 95 kDa

Preparation

This antibody was obtained by immunization of INSR recombinant protein (Accession Number: NM_001079817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111851(INSR antibody) at dilution of 1:500



IP Result of anti-INSR (IP:Catalog No:111851, 4ug; Detection:Catalog No:111851 1:500) with HeLa cells lysate 1600ug.