

TAF9B Polyclonal Antibody



Catalog Number: 162851

Product name

TAF9B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TAF9B

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human TAF9B (NP_057059.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TAF9B Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using TAF9B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.