

Mouse TALLA-1 / TSPAN7 (CD231) (His Tag) recombinant protein



Catalog Number: 502851

General Information

Gene Name Synonym

Cell surface glycoprotein A15; PE31; TALLA homolog; Transmembrane 4 superfamily member 2

Protein Construction

A DNA sequence encoding the mouse TSPAN7 (Q62283) (Arg113-Met213) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Arg 113

Molecular Mass

The recombinant mouse TSPAN7 comprises 112 amino acids and has a predicted molecular mass of 13 kDa. The apparent molecular mass of the protein is approximately 23-29 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse TALLA-1 / TSPAN7 Protein (His Tag) SDS-PAGE