

MRPS24 antibody



Catalog Number: 112851

Product name

MRPS24 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MRPS24 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human brain tissue. Observed molecular weight by Western-blot: 19kd

Preparation

This antibody was obtained by immunization of MRPS24 recombinant protein (Accession Number: NM_032014). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112851(MRPS24 antibody) at dilution of 1:500