

Anti-ANGPTL7 antibody



Catalog Number: 175851

Product name

Anti-ANGPTL7 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ANGPTL7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human angptl7. aa range:297-346

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of SH-SY5Y cells using Angptl7 Polyclonal Antibody diluted at 1: 2000